

The load level		Human weight
Light	●	< 60 kg
Average	●	60 ÷ 90 kg
Heavy	●	> 90 kg

FLEXIBLE POLYURETHANE FOAM – RECOMMENDED APPLICATIONS FOR FURNITURE								
	Pillows	Armrests	Back cushion	Backs	Seat cushion	Seats	Mattress cushion	Mattress
<b>STANDARD POLYURETHANE FOAM</b>								
VB2520		● ● ●	● ●		●		●	
VB2535				●		●		●
VB2540						●		●
VB3030				● ●	● ●	● ●	● ●	● ●
VB3035		● ● ●	● ●	● ●	● ●	● ●	● ●	● ●
VB3040						●		●
VB3050						●		●
VB3533			●	●	●	●	●	●
VB3540				●		●		●
<b>COMBUSTION MODIFIED STANDARD POLYURETHANE FOAM</b>								
VB2530A				●		●		●
CV2533				●		●		●
CV3037						●		●
CV3538					●	●	●	●
<b>HIGH RESILIENCE POLYURETHANE FOAM</b>								
HR2313	● ● ●	● ● ●	●					
HR2618	● ● ●	● ● ●	● ●					
HR3010	● ● ●	● ● ●						
HR3020		● ● ●		●	● ●	● ●	● ●	
HR3030				●	● ●	● ●	● ●	
HR3532			●	●	●	●	●	●
HR3538				● ●	● ●	● ●		● ●
HR3737				● ●	● ●	● ●		● ●
HR4030			●	●	●	●	●	●
HR4037				●		●		●
HR4250						●		●
HR5050						●		●
<b>COMBUSTION MODIFIED HIGH RESILIENCE POLYURETHANE FOAM</b>								
CH2715		● ● ●	● ●					
CH3015		● ● ●	● ●					
CH3025				● ●	●		●	
CH3031				●		●		●
CH3420		● ● ●	●	●	●	●	●	●
CH3525			●	●	●	●	●	●
CH3532					●	●		●
CH4038						●	●	●

RECOMMENDED USE FIELD			
Seasonal Furniture	Sponges of various colors	Shoulder padding	Noise isolation, packing material
VB1820, VB1830, VB1840, VB2220, VB2230, VB2240, CV2130	VB1830, VB1840, VB2102 – sponge foam VB2501 - for massage	VB1840, VB2240	VB1820

\* Experience has shown that high resilience polyurethane foam (HR and CH), produced in "Vita Baltic International" adequately chosen and used in accordance with the Recommendations, the basic properties and a sufficient level of comfort retain up to 10 years, standard polyurethane foam (VB and CV) up to 7 years.

\*\* These recommendations inappropriate for the foam, which is intended for intensive use, for example in sports halls, hospitals, kindergartens, etc.